

Department of Health Wadsworth Center

GENERAL ASSAY APPROVAL

Please submit all information as outlined below. Submit one hard copy of the entire package and one electronic copy (as a PDF file on a CD or flash drive) to:

US Postal Service: Clinical Laboratory Evaluation Program, Biggs Laboratory, Wadsworth Center, New York State Department of Health, Empire State Plaza, Albany, NY 12237; Attn: Assay Validation Review **UPS, FedEx, Courier:** Clinical Laboratory Evaluation Program, Biggs Laboratory, Wadsworth Center, New York State Department of Health, Dock J - P1 Level, Empire State Plaza, Albany, NY 12237; Attn: Assay Validation Review

Materials submitted, including related data packages, will not be returned to the laboratory. All materials are maintained under strict confidentiality. Materials are subject to New York State's Freedom of Information Law (commonly called FOIL). We suggest marking your documents as "proprietary" or "confidential". If so marked, laboratories will be given an opportunity to block information release.

SECTION 1: GENERAL INFORMATION	
Laboratory Name:	NYS PFI:
Assay (Test) Name*:	
Assay (Test) Type: 🗆 Laboratory Developed Assay (LD	A) □ Investigational Use Only (IUO)
□ Analyte Specific Reagent (ASI	R) \Box Research Use Only (RUO)
For RUOs only, provide the manufacturer/kit used:	
Target Population (if applicable):	
Methodology (e.g., ELISA, MIA, IFA):	
Analyte(s) included (West Nile Virus IgM, etc.):	
Validated Specimen Type(s) (e.g., serum, plasma, CSF	-):
Clinical Purpose:	

*If the test comprises a number of individual assays that are combined into a panel and interpreted as such you must clearly describe the composition and application of the panel. Complete packages for each individual assay must be submitted. In addition to validation of each individual assay, the combined panel must also be validated in its combination. For example, if the result consists of a risk/prediction factor, score or similar, it is this value that must be shown to be accurate, precise and reproducible, and meet all the other criteria described in the general requirements for assay validation.

All submissions must include a "Risk Attestation Form" found on the CLEP Test Approval webpage.

If this submission is a modification of an FDA or NYS-approved assay, see the "Assay Modification Checklist" on the <u>CLEP Test Approval webpage</u>.

If this submission is for the addition of an assay under an approved exemption, see the "Add Test Under Exemption Checklist" on the <u>CLEP Test Approval Webpage</u>.

SECTION 2: INSTRUCTIONS FOR SUBMITTING A FULL VALIDATION PACKAGE

The checklist below is a guide for items that must be included in the full validation submission packages. The information submitted must be organized as **numbered or uniquely named attachments**. If an item is not included, indicate the reason. Refer to the New York State General System Standards and any relevant Specialty Standards in preparing your submissions.

Section 2.1: Standard Operating Procedure and Controls

Procedure manuals must contain all required elements as described in the NYS General Systems Standards, Test Procedure Content Standard of Practice 1 (TPC S1): Test Procedure Content

File Name

File Name	
	Practitioner and patient educational materials that include a description of assay limitations and, where applicable, other information as may be necessary for informed consent of test subjects.
	Clinical indications for testing, including, where appropriate, the prevalence and description of the medical condition.
	Test subject preparation, specimen collection and handling, specimen rejection criteria, including a description of the mechanism to assure collection and transport requirements have been followed.
	A description of the assay, assay principle and clinical validity. For molecular tests, a description of the structure of the gene(s) to be tested, if applicable, must be provided.
	Complete and detailed procedures for performing the assay, including algorithms and flowcharts as necessary and any safety considerations.
	List of equipment / instrumentation essential to the assay.
	Reagents: source, preparation, storage stability and handling (amplification assays: include list of primers and sequences).
	Source and verification of standards / calibrators, quality control materials and the type, number, frequency and placement of the QC samples in an analytical run. Include QC evaluation/monitoring protocols.
	Calculation of results and interpretation (amplification assays: describe product size and method used to confirm product and result, where applicable).
	Assay interferences and limitations.
·	

File Name

A sample requisition form containing all the required elements in Test Request Standard of Practice 3 (TR S3): Test Request Form.
Examples of test reports containing all findings (e.g., positive, negative, indeterminate, inconclusive, etc.) with interpretive text, assay limitations and any disclaimers required by the federal government for tests utilizing Analyte-Specific Reagents (ASRs) and in compliance with Reporting Standard of Practice 2 (REP S2): Test Report Content .
If preliminary or presumptive positive results will be reported without confirmation, include examples of these reports containing the appropriate statements explaining the presumptive/preliminary nature of the results and recommendations for confirmation, if applicable.

Section 2.3: References

File Name

	A list of relevant literature references that describe the scientific basis and clinical validity of the assay. Provide electronic copies of primary references only.
	Applicable package inserts

Section 2.4: Validation

File Name	Validation Studies
	NARRATIVE SUMMARY of the validation studies performed with results and conclusions must be submitted and must be supported by providing the laboratory's validation protocols. The summary must address how analytical and clinical performance characteristics were established and describe any comparative methods and the source and number of specimens. Raw data must be provided using an appropriate number of samples across all representative specimen matrices and expected outcomes. Data should be summarized with <u>clearly labeled</u> tables, figures, and photographs.
	Analyte and specimen matrix stability
	Specimen transport conditions
	Storage time and temperature
	Accuracy
	Precision (reproducibility, both within (intra-) and between (inter-) runs)
	Reportable range, where applicable (calibration of quantitative tests).
	Analytical specificity, address potential cross-reactivity (for infectious disease testing) and any interferences (endogenous and exogenous)
	Clinical validity (sensitivity and specificity) establishment:
	Describe the protocols used to determine the clinical status of test subjects
	Describe the procedure used to blind the clinical status of specimens during testing
	Describe the procedures used to resolve discrepant or equivocal test results
	Present data used for the determination of clinical sensitivity, specificity and/or predictive values
	A description of studies performed to validate any data reduction and interpretation processes, including statistical or algorithmic calculations.

	A description of how reference intervals or assignment of cutoff values were determined, if applicable. For molecular tests, provide high quality original results of a sample of validation data.
File Name Representative Specimen Run	
	Provide actual instrument printouts, worksheets, or charts from one representative run, including all calibration and quality control materials.
	For those methods using gels, blots and/or electrophoresis, submit high quality photographs.
File Name Quality Control Data	
	Provide quantitative QC data and QC review format (e.g., Levey-Jennings charts) for 20-30 consecutive runs.

Section 2.5: Quality Assurance

File Name

Proficiency Testing: A plan for proficiency testing including criteria for passing, the number of
samples included in the proficiency panel, and the corrective actions that would be taken. The
proficiency testing plan must be compliant with Proficiency Testing Standards of Practice.
If there are no commercially available proficiency testing programs for the analyte(s), a detailed plan for an alternative to proficiency testing must be provided. This plan must be in compliance with Proficiency Testing Standard of Practice 3 (PT S3): Alternative to Proficiency Testing.
Systems agreement: If different instruments or platforms will be used to perform the assay, demonstrate the assay's consistency across these variables.
Quality Control Plan: Identify the critical steps in the test procedure and the quality control measures taken to control and monitor assay performance for consistent and reliable results.

I